

NZYTech qPCR Selection Guide



NZYTech's qPCR products were developed for highly accurate real-time PCR compatible with a vast range of thermocyclers.

Please use the following selection guide to choose the product that is most appropriate for your equipment:

PASSIVE REFERENCE DYE

EQUIPMENT *	High	Low	None
Agilent			
AriaMX MX3000P™, MX3005P™, MX4000P™		8	
Applied Biosystems™			
7000/7300/7700 7500/7500 FAST 7900/7900HT/7900HT FAST QuantStudio™ 6, 7, 12k Flex/ViiA7™ StepOne™/StepOne™ plus	8 8 9	8	
Bio-Rad®			
CFX96 [™] /CFX384 [™] Opticon [™] /Opticon [™] 2			8
Fluidigm [®]			
BioMark™		•	
Illumina®			
Eco™			8
Qiagen			
Rotor-Gene™ 3000 Rotor-Gene™ 6000 Rotor-Gene™ Q			8
Roche			
Lightcycler® 96 Lightcycler® 480 Lightcycler® Nano (*) For other equipment please contact NZYTech	(info@nzvtech.com)	X	8
() To other equipment pieuse contact NZT Tech	(wiley ceeieoi 11)		



www.stratech.co.uk/nzytech



info@stratech.co.uk +44 (0) 1638 782600



in @stratech-scientific-ltd

(i) @stratech.scientific







www.stratech.co.uk



